

November 21, 2023

## URGENT MEDICAL DEVICE RECALL NOTICE NURSE ASSIST Sterile 0.9% Saline, USP and Sterile Water, USP

R2331 1<sup>st</sup> Attempt

Nurse Assist has advised The Stevens Company of a Type II recall for the following products sold by The Stevens Company:

Product Name	<b>Stevens Catalog</b>	Manufacturer #	Expiry Date
	#		
WATER STERILE 120ml STERILE	372-6210	6210	
FOIL LID	372-0210	0210	
SOLUTION SALINE 0.9% NACL 120ml	372-6220	6220	
STER IRRIGATION FOIL LID	372-6220C	0220	
SOLUTION SALINE 0.9% NACL 100ml	372-6240	6240	
IRRSTER SCREW TOP BTL	372-0240	0240	
WATER STERILE 100ml IRRG	372-6250	6250	
w/SCREWTOP BOTTLE	372-0230	0230	
WATER STERILE 250ml IRRIGATION	372-6260	6260	Any product with an
SUCTIONING SCREW TOP BOTTLE	372-6260-X	0200	expiry date of:
SOLUTION SALINE 0.9% NACL 250ml	372-6270	6270	
STER IRRIGATION & SUCTIONING	372-6270-X	0270	November 1, 2023
SOLUTION SALINE 0.9% NACL 500ml	372-6280	6280	through
STER IRRIG SCREW CAP BOTTLE	372-0260	0280	September 18, 2025
SOLUTION SALINE 0.9% NACL 1000ml	372-6281	6281	
STER IRRIG SCREW CAP BOTTLE	372-0261	0281	
WATER STERILE 500ml IRRIG	372-6290	6290	
BOTTLESCREW CAP	372-0290	0290	
WATER STERILE 1000ml IRRIG	372-6291	6291	
BOTTLEBOTTLE SCREW CAP	372-0291	0291	
SOLUTION SALINE 0.9% NACL 7.1oz	272 6507	6507	
STER SPRAY WOUND FLUSH	372-6507	6507	

## Reason for Recall

Nurse Assist is conducting a recall for select products expiring within the range of November 1, 2023 through September 18, 2025. Any product with an expiry date within this range is recalled.

Nurse Assist has determined that the sterility assurance level (SAL) cannot be verified to the required degree of 10<sup>-6</sup>.

Product with an affected expiration date is to be destroyed in the field and no return is required. Credit will be issued for recalled product destroyed.

Products with an unverified sterility may be contaminated which could result in infection. The Stevens Company has received no reports of adverse events because of this issue.

Continued on next page ...

## FRM.PU.PR-01.2

Toll Free: 800-665-0368

## Action Required (Due January 31, 2024)

- 1. Immediately review your inventory for product with an impacted expiry date.
- 2. Record the lot numbers and expiry dates of all recalled products on the Acknowledgment Form accompanying this notice.
- 3. Destroy recalled products.
- 4. Notify all affected users, facilities or customers of the contents of this notice.
- 5. Complete the attached Acknowledgment and return by email to qc@safecross.com or fax to 416-665-6897.

Even if you have no inventory of this product, please complete the Acknowledgment and check the appropriate box. Your response ensures that this recall is carried out effectively.

Thank you, Steve Baker, CQIA Regulatory Affairs & Quality Assurance Manager

